

ACUPOINT TREASURE



Based on ancient Chinese acupressure, electric pulse, and magnetic therapy. Our Acupoint Treasure combines three therapies to give you a natural healing tool.

BENEFITS:

Magnetic Therapy

The human body is made up of electromagnetic vibrations.

Magnetic Therapy involves the use of magnetic fields to help activate the body's natural impulses to heal.

Chinese Acupressure (also known as Shiatsu)

Traditional Chinese Medicine is based on the concept that when the body is in balance, yin and yang, health is the result. Any imbalance and disruption to the flow of energy (known as qi) along meridian pathways can lead to disease.

Electric Pulse Therapy

This therapy involves the use of electric pulses to temporarily open the membrane barrier to allow molecules to enter the cell.

Acupoint Treasure uniquely combines all three of these methods, to provide you with a tool to naturally balance your body and improve your health.

UNIQUE FEATURES & BENEFITS:

Magnetic Therapy

Dr. Kyochi Nakagawa of Japan believes that because the earth's magnetic field has decreased by as much as 30 percent over the last millennium, magnetic therapy can provide what is missing. He suggests that this is why many of our modern ailments exist, he calls this "Magnetic Field Deficiency Syndrome"

Chinese Acupressure

There are over 2,000 acupressure points on the human body that can connect with the 12 main meridians and 8 secondary meridians.

Acupuncture originated in China more than 2,000 years ago, making it one of the oldest and most commonly used medical procedures in the world.

DID YOU KNOW?

In India it is said the builders of the pyramids were well acquainted with magnetic forces and used it in the preservation of mummies from natural decay. It is said Cleopatra used a small magnet on her forehead to preserve her beauty.

The use of magnets therapeutically has increased in the last few years, with professional athletes using them to relieve joint pain, and as a result, magnetic pain relief products are increasingly being sold in golf shops.

In the past two decades, acupuncture has grown in popularity in the United States. The report from a Consensus Development Conference on Acupuncture held at the National Institutes of Health (NIH) in 1997 stated that acupuncture is being "widely" practiced--by thousands of physicians, dentists, acupuncturists, and other practitioners--for relief or prevention of pain and for various other health conditions.

(Culliton PD. Current utilization of acupuncture by United States patients. Abstract presented at: National Institutes of Health Consensus Development Conference on Acupuncture; 1997)

All material in this document is provided for educational purpose only and may not be construed as medical advice or instruction. No action or inaction should be taken based solely on the contents of this information document. Instead, readers should consult appropriate health professionals on any matter relating to their health and well-being. All information and opinions provided in this document are believed to be accurate and sound, based on the best judgment available to the author, but the readers who fail to consult with the appropriate health authorities assume the risk of any injuries. The publisher is not responsible for errors or omissions.

KEY SCIENTIFIC STUDIES:

1. A study was conducted to determine whether subjects with post-polio syndrome and muscular or arthritic pain would benefit from using a permanent active magnet. Subjects using the magnets reported significantly less pain than those using a placebo.
(Vallbona, Carlos et al., Response of pain to static Magnetic fields in postpolio patients, a double blind Pilot study, Archives of Physical Medicine and Rehabilitation, 1997)
2. Magnet therapy reduced pain intensity in fibromyalgia patients, after sleeping on pads with active magnets.
(Journal of Alternative and Complementary Medicine, 2001)
3. Preliminary research suggests that patients with chronic asthma who receive acupressure may experience improved quality of life
(Maa SH, et al., Effect of acupuncture or acupressure on quality of life of patients with chronic obstructive asthma: a pilot study. Altern Med, 2003)
4. Following chemotherapy, patients who were treated with acupressure experienced less nausea and vomiting. The effect usually lasted for 6-8 hours, but they found that by placing an elasticized plastic stud over the point and pressing, could prolong the effect for 24 hours.
(Prolongation of the antiemetic action of P6 acupuncture by acupressure in patients having cancer chemotherapy. Dundee JW; Yang J, JR Soc Med, 1990)
5. Acupressure has shown effectiveness against motion sickness
(Bertolucci LE, et al. Efficacy of a portable acustimulation device in controlling seasickness, Aviat Space Environ Med 1995)

DIRECTIONS OF USE:

Usage: Remember three words: Press, Scan and Pulse.

- A. Press: Press the pinpoint of the Acupoint Treasure on the specific acupoint, using the strength, between 100-500g, according to the highest level the patient can tolerate.
- B. Scan: While pressing with strength, circling the acupoint with pinpoint of the Acupoint Treasure clockwise or counter clockwise (like massage), and at the same time scanning the acupoint with the magnetic bunch to change the still magnetism into active magnetism and change part of the bio-electric field of the body and by this way, to speed up the blood flow, making the healing process more significant.
- C. Pulse: Pressing the acupoint 1-3 minutes and pulsing the acupoint once or twice can clear and activate the channels and collaterals, and get effectiveness quickly. Before electrical pulse install the transmitting stick according to illustrated instructions, the stick must touch the skin, and keep the pinpoint 1 mm from the acupoint.
- D. Using two Acupuncture Treasure the same time on two different acupoints or the symmetrical acupoints of the same pinpoint (eg: temple, Renying acupoint, acupoints on the eyes etc) can get better effect.
- E. Moving the Acupoint Treasure vertically on the Meridian and collateral can also get good effect. For the diseases, which can cause pains, pressing directly on the Ahshih acupoint (aching point) e.g.: periarthritis of humeroscapularis, headache, waist pain, and nerve pain. For all ailments, please consult and follow the relative acupoints and press the specific acupoint with the Acupuncture of the Chinese acupoint with the Acupuncture Treasure.

Points to note:

- A. Read the directions carefully before using the Acupuncture Treasure.
- B. Prevent the Acupoint Treasure from shaking, moisture, high temperature and heat source. Keep Acupoint Treasure away from hand watches, TV, video recorders, which are susceptible to magnetic force.
- C. Do not disassemble without authorization; otherwise, the precision structure such as magnetic wave will be destroyed. The useful life of the electric pulse is generally 5000 times/piece.
- D. Patients with cardiac pacemaker should be cautious when using Acupuncture Treasure.

All material in this document is provided for educational purpose only and may not be construed as medical advice or instruction. No action or inaction should be taken based solely on the contents of this information document. Instead, readers should consult appropriate health professionals on any matter relating to their health and well-being. All information and opinions provided in this document are believed to be accurate and sound, based on the best judgment available to the author, but the readers who fail to consult with the appropriate health authorities assume the risk of any injuries. The publisher is not responsible for errors or omissions.